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Pedalo® Footboards moveable fulcrum:

- ✓ Stabilization of foot and leg axle
- ✓ Improvement of posture and movement coordination

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Pedalo® Pedasan:

- ✓ Sensitization of reflex zones
- ✓ Relaxation of deep musculature

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Further information on the issue of „A healthy foot“ on our website www.pedalo.de or in the brochure with the topic „A healthy foot“.

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Our feet are the fundamental base of our body and their significance should not be underestimated. It is generally known that our feet reflex our whole body. Our feet have a very complex structure with many bones, tendons and muscles that are created in a wonderful way.

Feet carry us for a lifetime and are still too often neglected. In its fine and highly complex construction, they are unique and do not require any improvement. In their construction they are not equal with anything. We just notice them when they call for help. From the beginning of our standing and running property the entire foot formation strengthens itself, so we always have it available. Walking barefoot is now mostly no longer possible anywhere. But it can be replaced by standing or walking on pebbles. In this way it can be seen how sensitive our feet are. Daily and regular exercises of this type comply with reflexology massage. It is a health care of high quality and also ensures the health of our entire body. It depends on how you use your foot. Deliberately and with feeling, that everyone can feel his feet again.


Recommended by Brigitte Wilhelmer, foot
reflexology practice, www.fussreflexologie.eu

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Using the Pedalo massage-mat is a treat for your feet, it improves your health and well-being. Due to the spherical curvature of the Pedalo massage-mat, all stimulus points of the foot soles are stimulated at the same time. Standing on the Pedalo massage-mat for only 3 minutes each day are enough not only to release your foot musculature and consequently the deep musculature of the whole body, but also to strengthen the cross arch of foot and to prevent foot suffering. Blood circulation and purification of the legs is stimulated and the cycle is activated. This creates a pleasant feeling of warmth. In sport it is for rapid recovery from stress. The massage mat is ideal to carry everywhere you want and it is also rapidly applicable e.g. away on travel, after training, in the shower cubicle.

Specification:
Natural rubber, ø 25 cm, 600 g

Important information:
This product is made of natural rubber. Slight discoloration of material are related to production process and therefore do not constitute a reason for complaint.

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First one, then the other foot, then both feet together. Do each of the three exercises 1 minute. Stand in all exercises with slightly bent knees, do not tense your toes and bring your ankle joints in a relaxed position.

Right foot:
Put your foot on the massage mat. The outside of your foot should show to the middle of the mat.

Left foot:
Repeat the exercise with your left foot.

Both feet:
Put both feet parallel on the massage mat. The insides of your foot should surround to the middle of the mat.

You can repeat the exercise, but just 1 minute per foot position! At the beginning pressure pains are normal. Always stay relaxed and let the naps act on your foot sole.

While you are sitting you can put your feet on the mat as long as it is comfortable for you. You can also sit on the mat. So your neck, shoulders, arms and hands can relax.

English